

## Antibodies Product Datasheet

|                               |                                                                                                                                                                      |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | NF-1C Monoclonal Antibody                                                                                                                                            |
| <b>Catalog Number:</b>        | N0717                                                                                                                                                                |
| <b>Reactivity:</b>            | Human, Mouse, Rat, Bovine, Chicken, Pig                                                                                                                              |
| <b>Applications:</b>          | WB                                                                                                                                                                   |
| <b>Dilutions:</b>             | Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.                                                                                                 |
| <b>Modification:</b>          | Unmodified/ Total                                                                                                                                                    |
| <b>Source:</b>                | Monoclonal Mouse                                                                                                                                                     |
| <b>Purification:</b>          | Affinity purification                                                                                                                                                |
| <b>Immunogen:</b>             | Purified recombinant human NF-1C protein fragments expressed in E.coli.                                                                                              |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                              |
| <b>Gene Name:</b>             | NFIC                                                                                                                                                                 |
| <b>Protein Name:</b>          | Nuclear factor 1 C-type                                                                                                                                              |
| <b>Other Name:</b>            | NFIC; NFI; Nuclear factor 1 C-type; NF1-C; Nuclear factor 1/C; CCAAT-box-binding transcription factor; CTF; Nuclear factor I/C; NF-I/C; NFI-C; TGGCA-binding protein |
| <b>Human Swiss Prot No.:</b>  | P08651                                                                                                                                                               |
| <b>Formulation:</b>           | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.                                         |
| <b>Specificity:</b>           | NF-1C Monoclonal Antibody detects endogenous levels of NF-1C protein.                                                                                                |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                         |