

Antibodies Product Datasheet

Product Name:	PA28γ Monoclonal Antibody
Catalog Number:	N0723
Reactivity:	Human, Mouse, Dog, Pig
Applications:	WB, IF/ICC
Dilutions:	Western Blot: 1/1000 - 1/2000. Immunofluorescence: 1/100 - 1/500. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human PA28γ protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	PSME3
Protein Name:	Proteasome activator complex subunit 3
Other Name:	PSME3; Proteasome activator complex subunit 3; 11S regulator complex subunit gamma; REG-gamma; Activator of multicatalytic protease subunit 3; Ki nuclear autoantigen; Proteasome activator 28 subunit gamma; PA28g; PA28gamma
Human Swiss Prot No.:	P61289
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	PA28γ Monoclonal Antibody detects endogenous levels of PA28γ protein.
Storage and Stability:	-20°C/1 year