

Antibodies Product Datasheet

Product Name:	TBL1XR1 Monoclonal Antibody
Catalog Number:	N0754
Reactivity:	Human, Mouse, Rat, Bovine, Chicken, Dog, Pig, Zebrafish
Applications:	WB, IHC-p, IF/ICC
Dilutions:	Western Blot: 1/1000 - 1/2000. Immunohistochemistry: 1/500 - 1/1000. Immunofluorescence: 1/100 - 1/500. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human TBL1XR1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	TBL1XR1
Protein Name:	F-box-like/WD repeat-containing protein TBL1XR1
Other Name:	TBL1XR1; IRA1; TBLR1; F-box-like/WD repeat-containing protein TBL1XR1; Nuclear receptor corepressor/HDAC3 complex subunit TBLR1; TBL1-related protein 1; Transducin beta-like 1X-related protein 1
Human Swiss Prot No.:	Q9BZK7
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	TBL1XR1 Monoclonal Antibody detects endogenous levels of TBL1XR1 protein.
Storage and Stability:	-20°C/1 year