

Antibodies Product Datasheet

Product Name:	TORC1 Monoclonal Antibody
Catalog Number:	N0761
Reactivity:	Human, Mouse, Rat, Pig
Applications:	WB, FCM, IF/ICC
Dilutions:	Western Blot: 1/1000 - 1/2000. Flow cytometry: 1/100 - 1/200. Immunofluorescence: 1/100 - 1/500. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human TORC1 (N-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	CRTC1
Protein Name:	CREB-regulated transcription coactivator 1
Other Name:	CRTC1; KIAA0616; MECT1; TORC1; WAMTP1; CREB-regulated transcription coactivator 1; Mucoepidermoid carcinoma translocated protein 1; Transducer of regulated cAMP response element-binding protein 1; TORC-1; Transducer of CREB protein 1
Human Swiss Prot No.:	Q6UUV9
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	TORC1 Monoclonal Antibody detects endogenous levels of TORC1 protein.
Storage and Stability:	-20°C/1 year