

Antibodies Product Datasheet

Product Name:	KDM5C/Jarid1C/SMCX mouse mAb
Catalog Number:	N0766
Reactivity:	Human
Applications:	WB, ICC
Dilutions:	wb 1:1000 icc 1:150
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human KDM5C / Jarid1C / SMCX protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	kdm5c
Protein Name:	
Other Name:	DXS1272E; Histone demethylase JARID1C; JARID1C; JmjC domain containing protein SMCX; Jumonji AT rich interactive domain 1C; Jumonji, AT rich interactive domain 1C (RBP2 like); Jumonji/ARID domain-containing protein 1C; KDM5C; KDM5C_HUMAN; Lysine (K) specific demethylase
Human Swiss Prot No.:	P41229
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KDM5C / Jarid1C / SMCX and does not cross-react with related proteins
Storage and Stability:	-20°C/1 year