

Antibodies Product Datasheet

Product Name:	KDM1/LSD1 mouse mAb
Catalog Number:	N0768
Reactivity:	Human, Monkey
Applications:	WB, ICC, IP
Dilutions:	wb 1:1000 icc 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human KDM1/LSD1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	kdm1a
Protein Name:	
Other Name:	Amine oxidase (flavin containing) domain 2; AOF2; BHC110; BRAF35 HDAC complex protein BHC110; BRAF35-HDAC complex protein BHC110; FAD binding protein BRAF35 HDAC complex, 110 kDa subunit; Flavin-containing amine oxidase domain-containing protein 2; KDM 1; KDM1; Kd
Human Swiss Prot No.:	O60341
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KDM1/LSD1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year