

Antibodies Product Datasheet

Product Name:	Dynactin 1(N-term) mouse mAb
Catalog Number:	N0775
Reactivity:	Human
Applications:	WB,IP
Dilutions:	wb 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant Dynactin 1 protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	dctn1
Protein Name:	
Other Name:	150 kDa dynein associated polypeptide;150 kDa dynein-associated polypeptide;DAP 150;DAP-150;DAP150;DCTN 1;DCTN1;DCTN1_HUMAN;DP 150;DP-150;DP150;Dynactin 1 (p150 Glued (Drosophila) homolog);dynactin 1 (p150 glued homolog Drosophila);Dynactin 1;Dynactin sub
Human Swiss Prot No.:	Q14203
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Dynactin 1(N-terminus),and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year

N0775 | Dynactin 1(N-term) mouse mAb

Tech@assaybiotech.com | Order@assaybiotech.com

Research use only not intended for therapeutic/diagnostic procedures

Last updated December 2022