

Antibodies Product Datasheet

Product Name:	Eg5 mouse mAb
Catalog Number:	N0776
Reactivity:	Human
Applications:	WB,IP,ICC
Dilutions:	wb 1:1000 icc 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Eg5 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	kif11
Protein Name:	
Other Name:	EG5;HKSP;KIF11;KIF11_HUMAN;kinesin family member 11;Kinesin like protein 1;Kinesin like spindle protein HKSP;Kinesin-like protein 1;Kinesin-like protein KIF11;Kinesin-like spindle protein HKSP;Kinesin-related motor protein Eg5;KNSL1;MCLMR;Thyroid receptor
Human Swiss Prot No.:	P52732
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Eg5 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year