

Antibodies Product Datasheet

Product Name:	DBC1 mouse mAb
Catalog Number:	N0777
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ICC, IP
Dilutions:	wb 1:500 icc 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human DBC1 protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	kiaa1967
Protein Name:	
Other Name:	DBCCR1; Deleted in bladder cancer protein 1; Deleted in bladder cancer protein 1 repressor; FAM5A; Protein FAM5A.
Human Swiss Prot No.:	Q8N163
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DBC1 and does not cross-react with related proteins
Storage and Stability:	-20°C/1 year