

Antibodies Product Datasheet

Product Name:	HDAC4(N-term) mouse mAb
Catalog Number:	N0778
Reactivity:	Human
Applications:	WB,IP
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human HDAC4 protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	hdac4
Protein Name:	
Other Name:	EC 3.5.1.98;HA6116;HD 4;HD4;HDAC 4;HDAC A;HDAC4;HDAC4_HUMAN;HDACA;Histone Deacetylase 4; Histone Deacetylase A;KIAA0288.
Human Swiss Prot No.:	P56524
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HDAC4 and does not cross-react with related proteins
Storage and Stability:	-20°C/1 year