

Antibodies Product Datasheet

Product Name:	KMT3C/SMYD2 mouse mAb
Catalog Number:	N0779
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human KMT3C / SMYD2 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	smyd2
Protein Name:	
Other Name:	Histone methyltransferase SMYD2; HSKM B; HSKM-B; HSKMB; KMT3C; Lysine N-methyltransferase 3C; MGC119305; N lysine methyltransferase SMYD2; N-lysine methyltransferase SMYD2; SET and MYND domain containing 2; SET and MYND domain containing protein 2; SET and MYND dom
Human Swiss Prot No.:	Q9NRG4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KMT3C / SMYD2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year