

Antibodies Product Datasheet

Product Name:	KLF4 mouse mAb
Catalog Number:	N0784
Reactivity:	
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant mouse KLF4 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	klf4
Protein Name:	
Other Name:	Endothelial Kruppel like zinc finger protein;Epithelial zinc finger protein EZF;EZF;GKLF;Gut-enriched krueppel-like factor;KLF;KLF4;KLF4_HUMAN;Krueppel-like factor 4;Kruppel like factor 4 (Epithelial zinc finger protein EZF) (Gut enriched Krueppel like fa
Human Swiss Prot No.:	O43474
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Transfected Only.
Storage and Stability:	-20°C/1 year