

Antibodies Product Datasheet

Product Name:	POU5F1/OCT4 mouse mAb
Catalog Number:	N0785
Reactivity:	Mouse
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant mouse POU5F1/OCT4 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	pou5f1
Protein Name:	
Other Name:	MGC22487;Oct 3;Oct 4;Oct-3;Oct-4;OCT3;Oct4;Octamer binding protein 3;Octamer binding protein 4;Octamer binding transcription factor 3;Octamer-binding protein 3;Octamer-binding protein 4;Octamer-binding transcription factor 3;OTF 3;OTF 4;OTF-3;OTF3;OTF4; P
Human Swiss Prot No.:	Q01860
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of POU5F1/OCT4 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year