

Antibodies Product Datasheet

Product Name:	Sox2 mouse mAb
Catalog Number:	N0786
Reactivity:	Mouse
Applications:	WB,FC,ICC
Dilutions:	wb 1:1000 icc 1:150
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant mouse Sox2 protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	sox2
Protein Name:	
Other Name:	ANOP3;cb236;Delta EF2a;icc;MCOPS3;MGC148683;MGC2413;RGD1565646;Sex determining region Y box 2;SOX 2;Sox2;SOX2_HUMAN;SRY (sex determining region Y) box 2;SRY box containing gene 2;SRY related HMG box 2;SRY related HMG box gene 2;SRY-box 2;Transcription fac
Human Swiss Prot No.:	P48431
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Sox2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year