

Antibodies Product Datasheet

Product Name:	COX IV mouse mAb
Catalog Number:	N0787
Reactivity:	Human, Mouse, Rat, Hamster, Goat, Monkey
Applications:	WB, FC, ICC, IP, IHC
Dilutions:	wb 1:1000 icc 1:150 fcm 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	A synthetic peptide corresponding to carboxyl terminal residues of human COX IV
Concentration:	1 mg/ml
Gene Name:	cox iv
Protein Name:	
Other Name:	AL024441; COX 4; COX IV 1; COX IV; COX IV-1; Cox4; COX41_HUMAN; Cox4a; COX4B; COX4I1; COX4I2; COX4L2; COXIV; Cytochrome c oxidase polypeptide IV; Cytochrome c oxidase subunit 4 isoform 1 mitochondrial; Cytochrome c oxidase subunit 4 isoform 1, mitochondrial; Cytochrome
Human Swiss Prot No.:	P13073
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of COX IV and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year