

Antibodies Product Datasheet

Product Name:	IL-8 mouse mAb
Catalog Number:	N0790
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	wb 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human IL8 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	il8
Protein Name:	
Other Name:	(Ala-IL-8)77;(Ser-IL-8)72;3 10C;310C;AMCF 1;AMCF1;b ENAP;Bbeta thromboglobulin like protein;Beta thromboglobulin like protein;C-X-C motif chemokine 8;CXC chemokine ligand 8;CXCL 8;CXCL8;Emoctakin;GCP 1;GCP-1;GCP/IL-8 protein I;GCP/IL-8 protein II;GCP/IL-8
Human Swiss Prot No.:	P10145
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects over-expressed and recombinant levels of IL8 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year