

Antibodies Product Datasheet

Product Name:	PPP1A mouse mAb
Catalog Number:	N0791
Reactivity:	Human
Applications:	WB,IHC,ICC
Dilutions:	wb 1:1000 icc 1:300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human PPP1A protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ppp1ca
Protein Name:	
Other Name:	Alpha isoform serine threonine protein phosphatase PP1alpha 1 catalytic subunit; Catalytic subunit;EC 3.1.3.16;MGC15877;MGC1674;PP 1A;PP-1A;PP1A;PP1A_HUMAN;PP1alpha;PP2C ALPHA;PP2CA;Ppp1ca;Protein Phosphatase 2C Alpha Isoform;Serine threonine protein phos
Human Swiss Prot No.:	P62136
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PPP1A and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year