

Antibodies Product Datasheet

Product Name:	Protein Phosphatase 4C mouse mAb
Catalog Number:	N0792
Reactivity:	Human
Applications:	WB,IHC,ICC
Dilutions:	wb 1:200 icc 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Protein Phosphatase 4C protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ppp4c
Protein Name:	
Other Name:	PP X;PP-X;PP4;PP4C;PP4C;PP4C_HUMAN;PPH3;PPP4;ppp4c;PPX;protein phosphatase 4 (formerly X), catalytic subunit;Protein phosphatase 4 catalytic subunit;Protein phosphatase X; Protein phosphatase X;protein phosphatase X, catalytic subunit;Serine/threonine pro
Human Swiss Prot No.:	P60510
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Protein Phosphatase 4C and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year