

Antibodies Product Datasheet

Product Name:	SMC1A(N-term) mouse mAb
Catalog Number:	N0796
Reactivity:	Human
Applications:	WB,IHC,FC
Dilutions:	wb 1:1000 fcm 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human SMC1A(N-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	smc1a
Protein Name:	
Other Name:	Chromosome segregation protein SmcB;DXS423E;KIAA0178;MGC138332;Sb1.8;Segregation of mitotic chromosomes 1;SMC protein 1A;SMC-1-alpha;SMC-1A;SMC1 (structural maintenance of chromosomes 1 yeast) like 1;SMC1;SMC1 structural maintenance of chromosomes 1 like
Human Swiss Prot No.:	Q14683
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SMC1A (N-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year