

Antibodies Product Datasheet

Product Name:	pyruvate dehydrogenase (lipoamide) α 1 mouse mAb
Catalog Number:	N0799
Reactivity:	Human, Mouse
Applications:	WB, ICC
Dilutions:	wb 1:1000 icc 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Pyruvate Dehydrogenase protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	pdha1
Protein Name:	
Other Name:	mitochondrial; OPA_HUMAN; PDH; PDHA; PDHA1; PDHCE1A; PDHE1 A type I; PDHE1-A type I; PHE1A; Pyruvate Dehydrogenase (lipoamide) α 1; Pyruvate dehydrogenase complex, E1 α polypeptide 1; Pyruvate Dehydrogenase E1 α ; Pyruvate dehydrogenase E1 component subu
Human Swiss Prot No.:	P08559
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of pyruvate dehydrogenase (lipoamide) α 1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year