

Antibodies Product Datasheet

Product Name:	KAP1/TIF1 β mouse mAb
Catalog Number:	N0800
Reactivity:	Human
Applications:	WB, ICC, IHC, IP
Dilutions:	wb 1:1000 icc 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human KAP1 / TIF1 beta protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	trim28
Protein Name:	
Other Name:	E3 SUMO protein ligase TRIM28; E3 SUMO-protein ligase TRIM28; FLJ29029; KAP 1; KAP-1; KRAB associated protein 1; KRAB interacting protein 1; KRAB-associated protein 1; KRAB-interacting protein 1; KRIP 1; KRIP-1; KRIP1; Nuclear corepressor KAP 1; Nuclear corepressor KA
Human Swiss Prot No.:	Q13263
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KAP1 / TIF1 beta and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year