

Antibodies Product Datasheet

Product Name:	Mannose Phosphate Isomerase mouse mAb
Catalog Number:	N0802
Reactivity:	Human,Rat
Applications:	WB,ICC
Dilutions:	wb 1:1000 icc 1:300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Mannose Phosphate Isomerase protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	mpi
Protein Name:	
Other Name:	PMI1;CDG1B;FLJ39201;Mannose 6 phosphate isomerase;Mannose-6-phosphate isomerase; MANNOSEPHOSPHATE ISOMERASE;MGC94106;MPI;MPI_HUMAN;Phosphohexomutase;phosphomannose isomerase 1;Phosphomannose isomerase;PMI;PMI1.
Human Swiss Prot No.:	P34949
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Mannose Phosphate Isomerase and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year