

Antibodies Product Datasheet

Product Name:	GATA6 mouse mAb
Catalog Number:	N0804
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human GATA6 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	gata6
Protein Name:	
Other Name:	Gata binding factor 6;Gata binding protein 6;GATA-binding factor 6;Gata6;GATA6_HUMAN; Transcription factor Gata 6;Transcription factor GATA-6.
Human Swiss Prot No.:	Q92908
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GATA6 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year