

Antibodies Product Datasheet

Product Name:	Hsp70(C-term) mouse mAb
Catalog Number:	N0806
Reactivity:	Human
Applications:	WB,IP
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Hsp70(C-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	hspa1a
Protein Name:	
Other Name:	DAQB 147D11.1 001;FLJ54303;FLJ54370;FLJ54392;FLJ54408;FLJ75127;Heat shock 70 kDa protein 1;Heat shock 70 kDa protein 1/2;Heat shock 70 kDa protein 1A/1B;heat shock 70kDa protein 1A;Heat shock 70kDa protein 1B;Heat shock induced protein;heat shock protein
Human Swiss Prot No.:	P08107
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Hsp70(C-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year