

Antibodies Product Datasheet

Product Name:	EpCAM mouse mAb
Catalog Number:	N0807
Reactivity:	Human
Applications:	WB,IP
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant EpCAM protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	cd326
Protein Name:	
Other Name:	17 1A; 323/A3; Adenocarcinoma associated antigen; Adenocarcinoma-associated antigen; Antigen identified by monoclonal antibody AUA1; AUA1; CD326; CD326 antigen; Cell surface glycoprotein Trop 1; Cell surface glycoprotein Trop 2; Cell surface glycoprotein
Human Swiss Prot No.:	P16422
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EpCAM and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year