

Antibodies Product Datasheet

Product Name:	TIF1 γ /TRIM33 mouse mAb
Catalog Number:	N0810
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human TIF1 gamma / TRIM33 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	trim33
Protein Name:	
Other Name:	8030451N04Rik;AI413936;cb1085;DKFZp586K1123;E3 ubiquitin-protein ligase TRIM33;EC 6.3.2.-;Ectodermin;Ectodermin homolog;FLJ11429;FLJ32925;id:ibd2175;MGC136680;mKIAA1113;OTTHUMP00000013662;OTTHUMP00000013663;Protein Rfg7;PTC7;Ret fused gene 7;RET-fused ge
Human Swiss Prot No.:	Q9UPN9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TIF1 gamma / TRIM33 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year