

Antibodies Product Datasheet

Product Name:	BiP/GRP78 mouse mAb
Catalog Number:	N0811
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, IP
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human BiP/GRP78 protein expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	hspa5
Protein Name:	
Other Name:	78 kDa glucose regulated protein; 78 kDa glucose-regulated protein; AL022860; AU019543; BiP; D2Wsu141e; D2Wsu17e; Endoplasmic reticulum lumenal Ca(2+)-binding protein grp78; Endoplasmic reticulum lumenal Ca2+ binding protein grp78; FLJ26106; Glucose Regulated Prote
Human Swiss Prot No.:	P11021
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BiP/GRP78 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year