

Antibodies Product Datasheet

Product Name:	CHD3(C-term) mouse mAb
Catalog Number:	N0813
Reactivity:	Human
Applications:	WB,IHC,FC
Dilutions:	wb 1:1000 1:500-1:1000 fcm 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CHD3 (C-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	chd3
Protein Name:	
Other Name:	226;2600010P09Rik;AF020312;ATP dependent helicase CHD3;ATP-dependent helicase Chd3;CHD-3;Chd3;CHD3_HUMAN;Chd7;Chromodomain-helicase-DNA-binding protein 3;hZFH;MGC40857;Mi 2 autoantigen 240 kDa protein;Mi 2a;Mi-2 autoantigen 240 kDa protein;Mi2 ALPHA;Mi2-a
Human Swiss Prot No.:	Q12873
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CHD3 (C-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year