

Antibodies Product Datasheet

Product Name:	HAUSP/USP7 mouse mAb
Catalog Number:	N0819
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human HAUSP / USP7 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	usp7
Protein Name:	
Other Name:	Deubiquitinating enzyme 7;HAUSP;herpes virus associated;Herpes virus associated ubiquitin specific protease;Herpesvirus-associated ubiquitin-specific protease;TEF 1;tef-1;TEF1;Ubiquitin carboxyl terminal hydrolase 7;Ubiquitin carboxyl-terminal hydrolase 7
Human Swiss Prot No.:	Q93009
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HAUSP / USP7 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year