

Antibodies Product Datasheet

Product Name:	Human Serum Albumin mouse mAb
Catalog Number:	N0820
Reactivity:	Human
Applications:	WB,ELISA,sELISA
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified Human serum albumin.
Concentration:	1 mg/ml
Gene Name:	alb
Protein Name:	
Other Name:	ALB;ALBU_HUMAN;Albumin (32 AA);Albumin (AA 34);Albumin;Analbuminemia;Bisalbuminemia;Cell growth inhibiting protein 42;DKFZp779N1935;Dysalbuminemic hyperthyroxinemia;Growth inhibiting protein 20;HSA;Hyperthyroxinemia dysalbuminemic;PRO0883;PRO0903;PRO1341;
Human Swiss Prot No.:	P02768
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Serum Albumin and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year