

Antibodies Product Datasheet

Product Name:	ATP-Citrate Lyase(C-term) mouse mAb
Catalog Number:	N0825
Reactivity:	Human,Mouse,Monkey
Applications:	WB,ICC,FC
Dilutions:	wb 1:1000 icc 1:150 fcm 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human ATP-Citrate Lyase protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	acly
Protein Name:	
Other Name:	ACL;Acly;ACLY_HUMAN;ATP citrate (pro-S) lyase;ATP citrate lyase;ATP citrate synthase; ATP-citrate (pro-S-)-lyase;ATP-citrate synthase;ATPcitrate synthase;ATPCL;Citrate cleavage enzyme;CLATP;EC 2.3.3.8;OTTHUMP00000164773.
Human Swiss Prot No.:	P53396
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ATP-Citrate Lyase and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year

N0825 | ATP-Citrate Lyase(C-term) mouse mAb

Tech@assaybiotech.com | Order@assaybiotech.com

Research use only not intended for therapeutic/diagnostic procedures

Last updated December 2022