

Antibodies Product Datasheet

Product Name:	DsbA-Tag mouse mAb
Catalog Number:	N0826
Reactivity:	E.coli
Applications:	WB-re
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant E.coli. DsbA protein fragments over expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	dsba
Protein Name:	
Other Name:	DSBA like thioredoxin domain containing protein.
Human Swiss Prot No.:	P0AEG4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DsbA and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year