

Antibodies Product Datasheet

Product Name:	FAS(C-term) mouse mAb
Catalog Number:	N0829
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human FAS (C-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	fas
Protein Name:	
Other Name:	ALPS 1A;ALPS1A;APO 1;Apo 1 antigen;APO 1 cell surface antigen;Apo-1 antigen;APO1;Apo1 antigen;APO1 cell surface antigen;Apoptosis antigen 1;Apoptosis mediating surface antigen FAS;Apoptosis-mediating surface antigen FAS;APT 1;APT1;CD 95;CD 95 antigen; CD9
Human Swiss Prot No.:	P25445
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FAS (C-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year