

Antibodies Product Datasheet

Product Name:	MEK1 mouse mAb
Catalog Number:	N0830
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant human MEK1 protein.
Concentration:	1 mg/ml
Gene Name:	MEK1/2
Protein Name:	
Other Name:	AA589381; Dual specificity mitogen-activated protein kinase kinase 1; Dual specificity mitogen-activated protein kinase kinase 2; EC 2.7.12.2; ERK activator kinase 1; ERK activator kinase 2; FLJ26075; MAP kinase kinase 1; MAP kinase kinase 2; MAP2K1; MAP2
Human Swiss Prot No.:	P36507,Q02750
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MEK1/2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year