

Antibodies Product Datasheet

Product Name:	PP2A α and β mouse mAb
Catalog Number:	N0831
Reactivity:	Human, Mouse, Rat
Applications:	WB, IP
Dilutions:	wb 1:2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human full length PP2A beta protein expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	ppp2cb
Protein Name:	
Other Name:	PP2A beta; PP2A-beta; PP2AB_HUMAN; PP2A beta; PP2CB; Ppp2cb; Protein phosphatase 2 (formerly 2A), catalytic subunit, beta isoform; Protein phosphatase 2 catalytic subunit beta isozyme; Protein phosphatase 2, catalytic subunit, beta isoform; Protein phosphatase 2A c
Human Swiss Prot No.:	P62714
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PP2A alpha and PP2A beta.
Storage and Stability:	-20°C/1 year