

Antibodies Product Datasheet

Product Name:	PP1C mouse mAb
Catalog Number:	N0832
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human PP1C protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ppp1cc
Protein Name:	
Other Name:	EC 3.1.3.16;PP 1G;PP-1G;PP1G;PP1G_HUMAN;PP1gamma;PPP 1G;PPP1CC;PPP1CC protein;PPP1G; Protein phosphatase 1 catalytic subunit gamma isoform;Protein phosphatase 1C catalytic subunit;Protein phosphatase 1C subunit;Protein phosphatase 2C gamma isoform; Serine
Human Swiss Prot No.:	P36873
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PP1C and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year