

Antibodies Product Datasheet

Product Name:	Angiotensin Converting Enzyme 1 mouse mAb
Catalog Number:	N0835
Reactivity:	Mouse
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Angiotensin Converting Enzyme 1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ace
Protein Name:	
Other Name:	ACE 1;ACE;ACE T;ACE_HUMAN;ACE1;Angiotensin converting enzyme somatic isoform;Angiotensin converting enzyme testis specific isoform;Angiotensin I converting enzyme;Angiotensin I converting enzyme 1;Angiotensin I converting enzyme peptidyl dipeptidase A 1;
Human Swiss Prot No.:	P12821
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Angiotensin Converting Enzyme 1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year