

Antibodies Product Datasheet

Product Name:	Ikaros(C-term) mouse mAb
Catalog Number:	N0843
Reactivity:	Human
Applications:	WB,IP
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Ikaros(C-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ikzf1
Protein Name:	
Other Name:	CLL associated antigen KW 6;DNA-binding protein Ikaros;hIk 1;hIk-1;Hs.54452;IK1;Ikaros (zinc finger protein);IKAROS;IKAROS family zinc finger 1 (Ikaros);Ikaros family zinc finger protein 1;Ikzf1;IKZF1_HUMAN;LYF1;Lymphoid transcription factor LyF-1;PRO0758
Human Swiss Prot No.:	Q13422
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Ikaros(C-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year