

Antibodies Product Datasheet

Product Name:	FOXO1(C-term) mouse mAb
Catalog Number:	N0844
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human FOXO1(C-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	foxo1
Protein Name:	
Other Name:	FKH 1;FKH1;FKH1;FKHR;FKHR;Forkhead (Drosophila) homolog 1 (rhabdomyosarcoma);Forkhead (Drosophila) homolog 1 (rhabdomyosarcoma);Forkhead box O1;Forkhead box protein O1; Forkhead box protein O1A;Forkhead in rhabdomyosarcoma;Forkhead, Drosophila, homolog of
Human Swiss Prot No.:	Q12778
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FOXO1(C-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year