

Antibodies Product Datasheet

Product Name:	Hexokinase 1 mouse mAb
Catalog Number:	N0847
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Hexokinase 1 protein fragments expressed in E.coli
Concentration:	1 mg/ml
Gene Name:	hk1
Protein Name:	
Other Name:	BB404130; Brain form hexokinase; dea; EC 2.7.1.1; Glycolytic enzyme; HEXOKIN; Hexokinase PI; Hexokinase type I; Hexokinase, tumor isozyme; Hexokinase-1; Hexokinase-A; HK I; HK1; HK1 tb; Hk1-s; HK1-ta; HK1-tb; HK1-tc; HKD; HKI; HMSNR; HXK1; HXK1_HUMAN; mHk
Human Swiss Prot No.:	P19367
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Hexokinase 1 and does not cross-react with Hexokinase 2 and other proteins.
Storage and Stability:	-20°C/1 year