

Antibodies Product Datasheet

Product Name:	IkB- α (N-term) mouse mAb
Catalog Number:	N0848
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human IkB-alpha(N-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	nfkbia
Protein Name:	
Other Name:	I kappa B alpha; I kappa B alpha; I(Kappa)B(alpha); I(Kappa)B(alpha); I-kappa-B-alpha; IkappaBalpHa; IkB-alpha; IKBA; IKBA; IKBA_HUMAN; IKBalpHa; MAD 3; MAD 3; MAD3; Major histocompatibility complex enhancer binding protein MAD3; Major histocompatibility complex enhanc
Human Swiss Prot No.:	P25963
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IkB-alpha(N-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year