

Antibodies Product Datasheet

Product Name:	Aurora Kinase A mouse mAb
Catalog Number:	N0849
Reactivity:	Human, Monkey
Applications:	WB, IF/ICC
Dilutions:	wb 1:500 icc 1:50
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Aurora Kinase A protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	aurka
Protein Name:	
Other Name:	AIK;ARK-1;ARK1;AURA;AURKA;Aurora 2;Aurora A;Aurora Family Kinase 1;aurora kinase A;Aurora Related Kinase 1;Aurora-related kinase 1;Aurora/IPL1 Like Kinase; Aurora/IPL1 Related Kinase 1;Aurora/IPL1-related kinase 1;AURORA2;AYK1;BRAK;Breast Tumor Amplified
Human Swiss Prot No.:	O14965
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Aurora Kinase A and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year