

Antibodies Product Datasheet

Product Name:	B Raf mouse mAb
Catalog Number:	N0850
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human B Raf protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	braf
Protein Name:	
Other Name:	FLJ95109; 94 kDa B raf protein; B raf 1; B Raf proto oncogene serine threonine protein kinase; B-Raf proto-oncogene serine/threonine-protein kinase (p94); BRAF 1; BRAF; BRAF_HUMAN; BRAF1; cRmil; MGC126806; MGC138284; Murine sarcoma viral (v-raf) oncogene homolog B1; M
Human Swiss Prot No.:	P15056
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of B Raf and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year