

Antibodies Product Datasheet

Product Name:	IKK β mouse mAb
Catalog Number:	N0851
Reactivity:	Human,Rat
Applications:	WB
Dilutions:	wb 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human IKK β protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	ikbkb
Protein Name:	
Other Name:	I kappa B kinase 2; I kappa B kinase beta; I-kappa-B kinase 2; I-kappa-B-kinase beta; I κ BKB; IKK 2; IKK B; IKK beta; IKK-B; IKK-beta; IKK2; IKKB; IKKB_HUMAN; IMD15; Inhibitor of kappa light chain gene enhancer in B cells; Inhibitor of kappa light polypept
Human Swiss Prot No.:	O14920
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IKK β and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year

N0851 | IKK β mouse mAb

Tech@assaybiotech.com | Order@assaybiotech.com

Research use only not intended for therapeutic/diagnostic procedures

Last updated December 2022