

Antibodies Product Datasheet

Product Name:	AMACR(C-term) mouse mAb
Catalog Number:	N0852
Reactivity:	Rat
Applications:	WB, ICC
Dilutions:	wb 1:1000 icc 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human AMACR(C-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	amacr
Protein Name:	
Other Name:	2 arylpropionyl CoA epimerase; 2-methylacyl CoA racemase; 2-methylacyl-CoA racemase; Alpha-methylacyl CoA racemase; Alpha-methylacyl Coenzyme A racemase; Alpha-methylacyl-CoA racemase deficiency, included; Alpha-methylacyl-CoA racemase; Amacr; AMACR deficiency, i
Human Swiss Prot No.:	Q9UHK6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AMACR(C-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year