

Antibodies Product Datasheet

Product Name:	MECT1/Torc1 mouse mAb
Catalog Number:	N0854
Reactivity:	Human
Applications:	WB,FC,ICC,IP
Dilutions:	wb 1:1000 icc 1:300 fcm 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human MECT1 / Torc1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	crtc1
Protein Name:	
Other Name:	KIAA0616;CREB regulated transcription coactivator 1;CREB-regulated transcription coactivator 1;CRTC1;CRTC1_HUMAN;FLJ14027;KIAA0616;MECT1;Mucoepidermoid carcinoma translocated 1;Mucoepidermoid carcinoma translocated protein 1;TORC-1;TORC1;Transducer of CR
Human Swiss Prot No.:	Q6UUV9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MECT1 / Torc1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year