

Antibodies Product Datasheet

Product Name:	CDK5(N-term) mouse mAb
Catalog Number:	N0856
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ICC
Dilutions:	wb 1:500 icc 1:150
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CDK5(N-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	cdk5
Protein Name:	
Other Name:	Cdk 5; Cdk5; CDK5_HUMAN; Cell division protein kinase 5; Crk6; Cyclin dependent kinase 5; Cyclin-dependent kinase 5; Protein kinase CDK5 splicing; PSSALRE; Serine threonine protein kinase PSSALRE; Serine/threonine-protein kinase PSSALRE; Tau protein kinase II catal
Human Swiss Prot No.:	Q00535
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CDK5(N-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year