

Antibodies Product Datasheet

Product Name:	CDC2/CDK1 mouse mAb
Catalog Number:	N0857
Reactivity:	Human
Applications:	WB, ICC
Dilutions:	wb 1:100 icc 1:50
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CDC2/CDK1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	cdk1
Protein Name:	
Other Name:	Cdc 2;Cdc2;CDC28A;CDK 1;CDK1;CDK1_HUMAN;CDKN1;CELL CYCLE CONTROLLER CDC2;Cell division control protein 2;Cell division control protein 2 homolog;Cell division cycle 2 G1 to S and G2 to M;Cell division protein kinase 1;Cell Division Cycle 2 Protein;Cyclin D
Human Swiss Prot No.:	P06493
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CDC2/CDK1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year