

Antibodies Product Datasheet

Product Name:	GRK2 mouse mAb
Catalog Number:	N0858
Reactivity:	Human
Applications:	WB
Dilutions:	wb 1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human GRK2 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	adrbk1
Protein Name:	
Other Name:	ADRBK 1;ADRBK1;Adrenergic beta receptor kinase 1;Adrenergic beta receptor kinase 1; ARBK1_HUMAN;BARK 1;BARK;BARK1;Beta adrenergic receptor kinase 1;Beta ARK 1;Beta ARK1; Beta-adrenergic receptor kinase 1;Beta-ARK-1;FLJ16718;G Protein Coupled Receptor Kina
Human Swiss Prot No.:	P25098
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GRK2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year