

Antibodies Product Datasheet

Product Name:	BTK mouse mAb
Catalog Number:	N0859
Reactivity:	Human
Applications:	WB,IP
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human BTK protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	btk
Protein Name:	
Other Name:	Agammaglobulinaemia tyrosine kinase;AGMX 1;AGMX1;AT;ATK;B cell progenitor kinase;B-cell progenitor kinase;BPK;Bruton agammaglobulinemia tyrosine kinase;Bruton tyrosine kinase;Bruton's Tyrosine Kinase;Btk;BTK_HUMAN;IMD 1;IMD1;MGC126261;MGC126262;OTTHUMP000
Human Swiss Prot No.:	Q06187
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BTK and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year